



Mount Protector Technical Bulletin & Instructions

About: The Mount Medix Mount Protector is a product specifically formulated for the taxidermy industry to clean, condition and prevent infestation in mounts. This product is a must for infestation prevention. The product is designed to be used minimally during semi-annual or annual cleanings. The formulation contains high-quality hair and hide conditioners to protect the mount, siliconized oils for easy dusting and no residue formation, plus it contains a small amount of the Mount Medix active insect ingredient to keep mounts from being damaged by moths and dermestids. The active will eliminate the mount as a food source for insects and their larvae. Mount Protector will restore the natural sheen to both hair and feathers, giving the mount a life-like appearance.

Mount Protector was created specifically for use by mount owners, regardless of familiarity with taxidermy, or taxidermy professionals.

Mount Protector is formulated to be easy to be a convenient, easy to use, pleasant smelling product to protect, clean and restore taxidermy mounts.

Instructions for use: Mount Protector is simply applied to all the mount surfaces, then wiped with a soft towel to remove excess. After cleaning and drying the mount will have been returned to its original condition, with renewed natural lustre and sheen, without looking wet or oily. Most importantly, **Mount Protector** will provide **protection** for insect infestation and sheds dust and dirt when used for routine cleaning. Mount Protector will leave the mount smelling fresh and clean.

Mount Protector can be used on most mount surfaces, including **hair, feather, antler, horns, hoof, paws and even habitat**. Mount Protector can provide years of worry-free protection against insect infestation.

Note: We do not recommend using Mount Protector on fish without testing its compatibility with the clear coats and finishes on the fish prior to use.

Applications:

Professional Use – Applying a coat or two of Mount Protector to the mounts prior to delivery will help keep your work looking good while in your client's care. This is accomplished by misting on a light coat, allow it to soak in for ten minutes and wipe with a soft cloth. Professional Taxidermists can also purchase the product wholesale, so they have it as a service to their clients.

Client Use – After semi-annual or annual cleaning with Mount Protector mist on a light coat, allow it to soak in for ten minutes and wipe with a soft cloth. If your taxidermist hasn't already made an initial application, let the Mount Protector thoroughly dry and apply a second light mist coat.

Special Application Note – Dedicated semi-annual or annual cleanings with Mount Protector will provide your taxidermy mounts indefinite protection from insect damage. That protection only lasts if you are dedicated in the use of Mount Protector.

Available sizes and usage:

The Mount Medix Mount Protector is available in 0,473-liter, 0,946-liter and 3,78-liter bottles Ready To Use sizes with an accompanying sprayer. Typically, a life-size Black Bear will require 1 – 1,5 liter for treatment. A typical shoulder mount Fallow Deer will require about 100-120 ml (1 - 1,2 dl). How many mounts it treats is heavily dependent on the size of the mount, thickness of hair and product absorption.

Shelf Life: Greater than 5 years

For additional information or questions please contact Mount Medix at info@mountmedix.eu